



CLEANING MANAGEMENT SYSTEM

#9 Multi-Purpose Disinfectant is designed to be used with the Cleanworks® Cleaning Management System. The system dispenses precisely diluted disinfectant/cleaner with the touch of a button...ready to work in spray, wipe-on or mop-on applications. Highly concentrated formula saves storage space and helps reduce disposal costs. Cleanworks® is simple, accurate and user friendly.

The cleanworks® Cleaning Management System is subject of U.S. Patents No. 5037003, 5083677, & 5228598

1.800.527.5722
www.cleanworks.net

9 MULTI-PURPOSE DISINFECTANT

Why it's a better cleaner

The Cleanworks® System dispenses it right...every time.

- Broad spectrum efficacy - Kills a wide range of organisms including HIV-1, MRSA, HBV and HCV.
- Multi-purpose detergent - Cleans soils from many hard surfaces
- Quaternary based - Hospital use disinfectant
- Water based, nonalkaline formula - Can be used on inanimate washable surfaces



How it cleans



Dispense precisely diluted cleaner into a clean, empty container. If solution is for reuse, properly label RTU container.



CLEAN DISINFECT :
Remove gross soil mechanically before cleaning begins.



Apply solution with wiper, mop or sponge, making sure to wet all surfaces thoroughly. Scrub. Discard solution when it becomes dirty and replace with fresh solution.



Allow surfaces to remain wet for 10 minutes. Remove excess liquid and let air dry.



Dispense 1,285 quarts or 321 gallons of ready-to-use disinfectant cleaner from the 1.25 gallon concentrate of Cleanworks® #9 Multi-Purpose Disinfectant.

Dispense 972 quarts* or 243 gallons** of ready-to-use disinfectant cleaner from the 2 liter concentrate case of Cleanworks® #9 Multi-Purpose Disinfectant.

Where it cleans

- Floors, Walls
- Mirrors, Glass
- Metal Surfaces
- Conductive Flooring
- Painted Surfaces
- Exterior Bowl Surfaces
- Empty Basins
- Non-Critical Medical Devices
- Hair/Nail Care Instruments
- Showers, Lavatory Fixtures
- Tanning Beds
- Other Hard, Nonporous Surfaces

Technical Data

Packaging (concentrate) 1.25 Gallon, 2/2 Liter

Concentrate Color Dark Blue

RTU Color Light Blue

RTU Fragrance Fresh

Foaming Low-Moderate

Free-Rinsing Yes

% Phosphates None

pH (concentrate) 7.2 - 8.2

EPA Registered Yes

Cleanworks® 2 RTU Yield per case

* Low-Flow 972 Quarts

** High-Flow 243 Gallons

Efficacy

Pseudomonas aeruginosa, *Staphylococcus aureus*, *Salmonella choleraesuis*, *Escherichia coli*, *Staphylococcus aureus* (MRSA), *Staphylococcus aureus* (VISA), *Streptococcus pyogenes*, *Enterococcus faecalis* (VRE), *Canine distemper*, *Feline leukemia*, *Hepatitis B* (HBV), *Hepatitis C* (HCV), *Herpes simplex type 1 & 2*, *HIV-1* (AIDS virus), *Pseudorabies*, *SARS Associated Coronavirus*, and many more.